


How To Handle UV-lamp with Chip DSC



Linseis Messgeräte GmbH

Gerlach

Stand: 18.03.2020

Index

1. General information.....	3
2. Prepare the Hardware	3
2.1 Requirements.....	3
2.2 assemble UV-Lamp.....	4
2.3 connect the UV-Lamp.....	4
3. Prepare and setup the measurement.....	5
3.1 setup the UV-option in software.....	5
3.2 setup measurement conditions in Acquisition software.....	5
4. Evaluate your measurement	6

1. General information

This manual is a short description for handling UV-lamp function for Chip DSC. The presented applications are just a few to show some possibilities. For more Information, read the other available instructions about software or specific manuals for the Chip DSC.

Important safety instruction:

Make sure all safety accessories and functions work. Wear protective glasses for UV-light all the time you work with Chip-DSC UV-lamp option! UV-radiation can harm eyes and skin.


2. Prepare the Hardware

2.1 Requirements

- Check if UV-Lamp and UV-Lamp power supply are on hand
- Check if USB Dongle with unlocked UV-Function is on hand
- Check if Chip-DSC UV-lamp measurement cell is on hand
- UV sensitive curing adhesive or other sample of interest is available


2.2 assemble UV-Lamp

- Install the UV-lamp in the fitting holder and fix it with screws
- Insert the UV-lamp holder in the UV-Lamp measurement cell
- Connect the safety switch cable with UV-lamp measurement cell
- Connect the UV-lamp with the UV-lamp power supply
- Connect the safety switch cable with the shown pins of UV-lamp power supply



2.3 connect the UV-Lamp

- Connect the RS232-USB connector with UV-lamp power supply and your PC
- Connect UV-lamp power supply with power net



- Turn on UV-lamp power supply


3. Prepare and setup the measurement

3.1 setup the UV-option in software

- Start the administration software and create a new instrument
- Choose “Uvitron skybeam” and press Ok
- Check the ComPort of your UV-lamp power supply (device manager → com ports) and insert this in device settings (administration software) of your UV-Lamp
- Insert a name for your UV-Lamp in device settings (administration Software) of your UV-lamp

3.2 setup measurement conditions in Acquisition software

- Start the acquisition software
- Setup your general measurement settings like name, weight and temperature profile
- Go to “Gas” Tab. If gas tab is not available go to the administrator software and create UV light option.



- Here you can setup your UV-light profile by clicking “+” and check and compare this with your temperature profile
- By double click on a single part of your profile the following window will appear where you can setup the settings to your needs. The UV lamp intensity can be adjusted in terms of % in the “nominal flow” field.

● Ramp segment

Nominal flow: 10,00 l

Rise time: 1,0000 min

Dwell: 1,000 min

○ Formula segment

Formula:

Dwell: 1,00 min

Flow change

Change nominal value of this segment to: 0,00 l

if the channel: Temperature

has a current value: higher

than: 0 °C

OK Cancel

- If you finished your setup start your measurement

4. Evaluate your measurement

- If your measurement finished start evaluation software and load your measurement
- In channel selection choose UV light and all other channels you are interested in
- With data explorer you can check and mark the power at a certain state in your measurement

