

SPECIFICATIONS

Optical Dilatometer complete with soft- and hardware, incl. reference samples and accessories

Data acquisition software for controlling of the Optical Dilatometer with several Parameter adjustments incl. temperature programming, gas control etc.

Optical Dilatometer L74/V/2000

Temperature range:	RT up to 2000°C
Sample length:	up to max. 20mm
Resolution:	up to 1µm
Atmosphere:	Inert, red., vac.
Vacuum:	up to 10E-2mbar
Reference sample:	reference sample with table included