

## SPECIFICATIONS

Optical Dilatometer complete with soft- and hardware, incl. reference samples and accessories

Data acquisition software for controlling of the Optical Dilatometer with several Parameter adjustments incl. temperature programming, gas control etc.

### Optical Dilatometer L74/HS/1700

<b>Temperature range:</b>	RT up to 1700°C
<b>Sample length:</b>	up to max. 20mm
<b>Resolution:</b>	up to 1 µm
<b>Vacuum:</b>	up to 10E-4mbar
<b>Referenz sample:</b>	referenz sample with table included