

## SPECIFICATIONS

Optical Dilatometer complete with soft- and hardware, incl. reference samples and accessories

Data acquisition software for controlling of the Optical Dilatometer with several Parameter adjustments incl. temperature programming, gas control etc.

### Optical Dilatometer L74/HS/500

<b>Temperature range:</b>	-100 up to 500°C
<b>Sample diameter:</b>	Max. 20mm
<b>Resolution:</b>	up to 1 µm
<b>Vacuum:</b>	up to 10E-4mbar
<b>Calibration standards:</b>	Calibration standard with table included
<b>Sample length:</b>	Up to max. 20mm