

SPECIFICATIONS

Optical Dilatometer complete with soft- and hardware, incl. reference samples and accessories

Data acquisition software for controlling of the Optical Dilatometer with several Parameter adjustments incl. temperature programming, gas control etc.

Optical Dilatometer L74/HS/600LT

| | |
|-------------------------------|--|
| Temperature range: | -150°C up to 600°C |
| Sample length: | Up to max. 20mm |
| Sample diameter: | Max. 7mm |
| Resolution: | up to 1 µm |
| Vacuum: | up to 10E-4mbar |
| Calibration standards: | Calibration standard with table included |
| Sample length: | Up to max. 20mm |